



Anti-bedsore mattress with variable pressure pump FDA P06 (Reconditioned)

Used, very good technical condition,
Production P.R.C.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Pressure pneumatic mattress: 70 - 120 mmHg,
Efficient quiet pump: 14 l / min,
Cycle time: 6/9/12/25 min,
Number of chambers: 20,
Usable area around: 215 x 105 cm,
Has the function:
ventilation
pillows,
CPR 'Quick Deflate' for emergency procedures,
Maximum patient weight: 210 kg,
Dimensions of the pump: 29 (L) x 12 (S) x 10 (H) cm,
Pump weight: 4 kg,
Mattress weight: 10 kg,
Included:
Bag for transport,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Installments, Leasing, Loan

