



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Bariatric bed MMO 8000 (New)

New
Polish production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Maximum patient weight: 250 kg,
Safe working load: 340 kg,
Lair filled with panels made of HPL board with a thickness of 12 mm,
Bed width: 1000/1150/1200 mm,
Width of the bed with side barriers: 1150/1250/1350 mm,
Bed length: base - 2000 mm, with extension up to 2180 mm,
Overall length: 2140 mm - 2320 mm,
Side barriers 480 mm high, with four crossbars, folded with one hand,
To ensure a better weight distribution, a set of eight double wheels with a diameter of 100 mm. From the top of the upper wheel, directional, central locking of all wheels activated from the bottom of the top,
Automatic locking of the central brake in the lowest position,
Electrical height control, backrest and thigh segment inclination angle, electrical CPR function. Implemented with the patient's wired driver and the control panel for the carer,
Height adjustment in the range: 210 - 830 mm,
Trendelenburg and anti-Trendelenburg position: 14 °,
Adjustment of inclination of the backrest segment: 70 °, 30 ° thigh segment, 20 ° leg segment,
When adjusting the backrest segment, automatic pause at 30 ° and 45 ° angle,
Electric and manual CPR, color-coded levers located on both sides of the bed,
LED night light, located under the bed,
Shelf for bedding,
Battery powered,
Drip stand,
Handles for fixing patient restraints,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Installments, Leasing, Loan

