



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Suction pump for M (ASPINA) (portable suction unit) (New)

- New
- Slovak production,
- **Technical data:**
 - Rated voltage / frequency:
 - 230V / 50 (60) Hz
 - 110V / 60Hz
 - Max. Current: 3.4 (3.8) A, 7.3A
 - Maximum flow: 1100 lit.min-1
 - Maximum negative pressure: 12 kPa
 - Loudness: ≤48 dB [A]
 - Weight: 36 kg
 - Dimensions (W x D x H) *: 565x350x860 mm
- Has a valid Technical Passport issued,
- 12 months warranty,
- Financing option (Poland only): Installments, Leases, Loans

