



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **SUCTION PUMP FOR EKOM 2.1-10 (stationary mammal in a box) (New)**

New

Slovak production,

Mammal in a silencing housing size DK50.10S,

Mammal comes standard with a control module that starts the device after the 24V AC / DC signal flow from the unit. The box effectively reduces the noise of the mammal.

Technical data:

Rated voltage / frequency:

230 V / 50 (60) Hz

110V / 60 Hz

Max. Current: 3.3 (3.7) A / 7.2 A

Maximum flow: 1100 liters min<sup>-1</sup>

Maximum vacuum: 12 kPa

Control voltage: 24V / 0.1A DC / AC

Operating mode: S1-100%

Volume without the box: ≤58 dB [A]

Volume with a box: ≤46 dB [A]

Weight without box: 15 kg

Weight with case: 29 kg

Dimensions without box (W x H x W): 330x430x350

Dimensions with the box (sz x g x w) \* 420x525x620

Has a current Technical Passport issued,

12 months warranty,

Possibility of financing (only Poland): Installments, Leasing, Loan

