



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BORADSOUND ES6E123 head endovaginal (Reconditioned)

Used, perfect technical condition,

Taiwanese production (2014),

Probe name: EC123,

Cooperation with Esaote ultrasound machines such as: Picus / DU3 / 4 / Caris Plus / My Lab 15/20/25/30,

Type: microconvex broadband, multi-frequency,

Operating frequency range: 4 - 8 MHz,

Application: endovaginalne,

Has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (only Poland): Installments, Leasing, Loan

