



## Endoscopic camera OLYMPUS OTV-SX2 (Reconditioned)

- Used, very good technical condition,
- Japanese production,
- Power supply: 230 V,
- Frequency: 50-60 Hz,
- Power: 50 VA,
- TV type: PAL,
- Adjusting the color saturation: red and blue,
- WHITE BALANCE function (White balance),
- Aperture adjustment,
- Video enhancement,
- **Out-put:**
  - Video OUT,
  - 2 x Y/C,
  - RGB,
  - H.COPY,
  - SYSTEM,
  - VTR REMONT,
- **The kit includes:**
  - Image processor OTV-SX2C,
  - OLYMPUS OTV-SX2 Camera Head adapted for rigid endoscopes, placed in a portable case,
- Dimensions: 29.5 x 38 x 9.5 cm,
- Weight: 4 kg,
- It has exoised the current Technical Passport,
- Warranty:
  - 6 months for national field - Poland,
  - 1 month for international field,
- The possibility of financing (Polish only): Installments, lease, loan.



