



Endoscopic camera PENTAX ENDO-VISION 3000 (Reconditioned)

- Used, very good technical condition,
- German production,
- Power supply: 230 V,
- Frequency: 50-60 Hz,
- Power: 19 W,
- Color system: PAL / NTSC
- Microswitches for adjusting brightness, gain and white balance,
- **Has sockets:**
 - 2 x S- VHS,
 - Out-put Audio,
 - 2 x BNC,
 - RGB,
- **Includes IPX7 Camera Head:**
 - Image format type PAL,
- Cat 02,
- It has exoised the current Technical Passport,
- Warranty:
 - 6 months for national field - Poland,
 - 1 month for international field,
- The possibility of financing (Polish only): Installments, lease, loan.

