



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video gastroscope Olympus GIF-2T160 (Reconditioned)

- Used,
- Japanese production,
- 3.7 mm duct for suction
- 2.8 mm duct for Endo-Therapy accessories
- 12.6 mm Outer diameter of the probe.
- 120° field of view.
- Has an auxiliary water hole to ensure continuous visibility
- Compatible with CV-180, CV-160, CV-140 and CV-100 video processor
- Diameter: 12.6mm
- Working length: 103 cm
- angulation (210 ° up, down and 90° 100° right / left)
- Has a valid Technical Passport issued,
- 6-month warranty,
- Financing option (Poland only): Installments, Leases, Loans,

VIEW IMAGES OF THE DEVICE

