



Video colonoscope OLYMPUS Evis Exera PCF-160AL (Reconditioned)

- Used, perfect condition,
- Japanese production,
- View direction: forward view,
- Outer diameter: 12 mm,
- The outer diameter of the tip: 11.5 mm,
- The diameter of the working channel: 3.2 mm,
- Working length: 1680 mm,
- Overall length: 2000 mm,
- Depth of field: 3-100 mm,
- Max. The angle of bending the tip of the endoscope up and down: 180° / 180°
- Max. The angle of bending the tip to the right and left: 160° / 160°,
- Maximum angle of field of view: 140°,
- Fully compatible with CV-180/160/140 processors,
- In set:
 - Original suitcase with dimensions: 59 x 43 x 16 cm,
 - Valves,
- User manual in English,
- Kat 01,
- Has an up-to-date review and Technical Passport,
- 6-month warranty,
- Financing option (Poland only): Installments, Leases, Loans,

VIEW IMAGES OF THE DEVICE

