



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Video duodenoscope OLYMPUS Evis Exera TJF-160VR (Reconditioned)**

- Used, perfect technical condition,
- Japanese production,
- View direction: 5 ° side view,
- Outer diameter:  $\varnothing$  11.3 mm,
- The diameter of the working channel: 4.2 mm,
- Working length: 1240 mm,
- Total length: 1550 mm,
- Depth of field: 5-60 mm,
- Air flow rate: 25 cm<sup>3</sup> / s,
- Max. The angle of bending the tip of the endoscope up and down: 120 ° / 90 °
- Max. The angle of bending the tip left and right: 110 ° / 90 °,
- Maximum angle of field of view: 100 °,
- User manual in English (PDF),
- Has a valid Technical Passport issued,
- 6-month warranty,
- Financing option (Poland only): Installments, Leases, Loans,

### **VIEW IMAGES OF THE DEVICE**

