



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video-ureterscopes OLYMPUS WA50200A (Ex-Display)

- New, exhibition,
- Japanese production,
- Direction of view: 7°,
- Tip lenght: 430 mm,
- Fiber lenght: 3000 mm,
- Canal: 4,2,
- Atraumatic tip,
- Extremely narrow introductory pipe,
- Image size increased by 25%,
- Port with two irrigation tips,
- Compatible with processor OTV S7,
- **Included:**
 - Tray OLYMPUS TRAY WA59383A,
- Dimensions: 59,5 x 29,5 x 5,5 cm,
- It has exoised the current Technical Passport,
- Warranty:
 - 6 months for national field - Poland,
 - 1 month for international field,
- The possibility of financing (Polish only): Installments, lease, loan.

