



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Cholescope OLYMPUS CHF-P20 (Reconditioned)

- Used, very good technical condition,
- Japanese production,
- The camera has one blind spot,
- View direction: forward view,
- Working length: 38 cm,
- Outside diameter: 4.9 mm,
- The diameter of the working channel: 2.2 mm,
- Max. The angle of bending the tip of the endoscope up and down:
 - Up: 160 °,
 - Down: 130 °,
- Maximum angle of field of view: 120°,
- Includes original case with dimensions: 59 x 43 x 16 cm,
- Has an up-to-date review and Technical Passport,
- 6-month warranty,
- Financing option (Poland only): Installments, Leases, Loans,

