



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video bronchoscope OLYMPUS BF-1T160 (Reconditioned)

- Used, perfect condition,
- Japanese production,
- View direction: forward view,
- Outer diameter: 6 mm,
- The outer diameter of the tip: 6 mm,
- The diameter of the working channel 2.8 mm,
- Working length: 600 mm,
- Total length: 870 mm,
- Depth of field: 3-100 mm,
- Max. The angle of bending the tip of the endoscope up and down: 180o / 130o
- Maximum angle of field of view: 120o,
- Fully compatible with CV-100/140/160 systems
- **In set:**
 - Original suitcase with dimensions: 59 x 43 x 16 cm,
- User manual in English,
- Has an up-to-date review and Technical Passport,
- Warranty: 6 months,
- Financing option (Poland only): Installments, Leases, Loans,

VIEW IMAGES OF THE DEVICE

