



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Light Source SOLOS ELS-2 (Reconditioned)

- Used, very good technical condition
- American production,
- Adapted to rigid endoscopes,
- Xenon light,
- 230 V power supply,
- 50-60 Hz frequency,
- 175 W xenon bulb,
- Digital display
- Has functions:
 - Smooth light intensity regulation,
 - Manual / auto,
 - 75 Ω,
 - Video In / Out,
- Equipped with carrying handles,
- Dimensions: 41 x 41x 20 cm.
- Has a valid Technical Passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 1 month on the international market,
- Financing option (Poland only): Installments, Leases, Loans.

