



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LCD Medical Monitor STORZ 200902 37 (Reconditioned)

- Used, very good technical condition,
- German production,
- Power supply: 230 V, through the power supply DC 12V,
- Frequency: 50/60 Hz,
- Power: 48 W,
- Image size: 19',
- Video input connectors: DVI-D,
- Dimensions: 447 x 387 x 85 mm,
- Weight: 8,5 kg,
- Cat 01,
- It has exoised the current Technical Passport,
- Warranty:
 - 6 months for national field - Poland,
 - 1 month for international field,
- The possibility of financing (Polish only): Installments, lease, loan.

