



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE CDA19 medical monitor (Reconditioned)

- Used, very good technical condition,
- Taiwanese production,
- Model Number: USE1913A,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power consumption: 100 W,
- 19 "color LCD display,
- Resolution: 1280 x 1024 pixels,
- Aspect ratio: 4: 3,
- Refresh: 60 Hz,
- Anti-reflective coating that minimizes reflection,
- Height adjustment in the range of 47 - 52 cm,
- Monitor angle adjustment (tripod): up / down,
- Removable cover covering the entrance ports,
- **Screen Settings:**
 - Brightness,
 - Contrast,
 - Colour,
 - Temperature,
- **I / O connectors:**
 - DVI
 - VGA,
- Includes power supply with power cord,
- Dimensions: 43 x 47.5 x 25 cm,
- Weight: about 4.5 kg,
- Has a valid Technical Passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 1 month on the international market,
- Financing option (Poland only): Installments, Leases, Loans

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

