



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynaecological chair SELENE M (New)

New,
Made in Poland,
It has a stable steel construction and fixed height,
The ABS base has a folding metal step for easy entry and exit,
Backrest and seat are adjustable with gas springs,
As a standard the armchair is equipped with a set of knee pads adjustable in all planes, a stainless steel bowl and a set of euroshins mounted on the sides of the seat,
There is a free option to order upholstery in the colour of your choice, except for "Presitige" colours which are available at an extra charge,
For an additional charge:
Paper hanger,
Removable footrest,
Hand grips (armrests),
IV hanger,
Injection support,
Technical parameters:
Length without footrest 139 cm, with footrest 190 cm,
Width 65 cm,
Height 85 cm,
Trendelenburg position 0° - 20°,
Anti-Trendelenburg position -15° - 0°,
Adjustment angle of backrest 0° - 75°,
Step height 25 cm,
Weight 70 kg,
Maximum load 160 kg,
It has a valid Technical Passport,
Warranty on structure 24 months,
Warranty on hydraulic accessories 12 months,
Possible financing (Poland only): Installments, Leasing, Loan.
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

