



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRAISTECH Goepel footrests (knee socks) (New)

New

Polish production,

They are used to support the lower limbs in a lying or sitting position,

They are equipped with polyurethane mattresses with dimensions: 17.5 x 27.5 x 3 cm,

The structure is made of stainless steel,

Footrests mounted on a ball joint allow adjustment in many planes,

The kit includes two footrests with arms and two mounting bones for the rail with adjustable height and angle,

Has a valid Technical Passport issued,

12 months warranty,

Financing option (Poland only): Installments, Leases, Loans,

