



OLYMPUS LF-GP laryngofiberscope (Reconditioned)

Used, very good technical condition,
Japanese production,
View direction: forward view,
Outside diameter: 4,1 mm,
The outer diameter of the tip: 3.8 mm,
Diameter of the working (suction) channel: 1.5 mm,
Working length: 600 mm,
Overall length: 855 mm,
Depth of field: 4 - 50 mm,
Maximum bending angle of the endoscope tip up and down: 120 ° / 120 °,
Maximum angle of view: 90 °,
Includes a suitcase with the dimensions: 33 x 9 x 40 cm,
Cat 02,
Has a valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month on the international market,
Financing option (Poland only): Installments, Leases, Loans,

