



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HEYINOVO WF-2423DG hot biopsy forceps (Ex-Display)

New, exhibition,
Reusable,
Autoclavable: 134 ° C,
For working channel with diameter: 2.8 mm,
Pipe diameter: 2.4 mm,
Working length: 2300 mm,
Expiry date: 03-2019,

