



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SURTRON TOUCH 200 surgical diathermy (New)

New
Italian production,
It is intended for small and medium-sized surgical procedures
It is suitable for precise monopolar and bipolar operations as well as microsurgery,
It enables deep personalization of settings and remembering for any number of operators.
The ability to switch the monopolar output functions as well as the power supply in the electrode holder allows you to perform surgery without distracting the surgeon from the place of surgery.
The level of contact safety between the separable passive electrode and the patient is controlled.
The neutral electrode connection is constantly monitored,
7 "touch screen that allows you to control device parameters.
Allows for minimally invasive surgery.
The electronic reading of the delivered power together with the monitoring of operational functions by the microcontroller guarantees complete reliability of the working conditions.
Thanks to its capabilities, it allows:
CUTTING,
Mixing-coagulation MIXED CUTTING,
Superficial coagulation surface,
Deep coagulation without dying of the tissue SOFT COAGULATION,
BIPOLAR coagulation,
BIPOLAR cut,
Application:
Outpatient, pediatric, vascular, plastic surgery,
Endoscopy,
Dermatology,
Otolaryngology,
Urology,
Gynecology,
Veterinary medicine,
First aid,
Gastroenterology,
Pneumologia,
Neurosurgery
Capabilities:
Can be used for operations with minimal invasiveness
Possibility to use bipolar pliers
The possibility of using a separate passive electrode
Planning of working conditions
Remembering the last settings
Switching on the monopolar electrode via a foot switch or located in the handle
Digital regulation and output power indication
Adjusting the intensity of information sounds
High performance in a small size
Activation in a handgrip or foot switch
Bipolar coagulation for safer treatments
Equipment:
Handle with switches,
Electrode set (10 pcs),
Passive electrode cable,
Steel passive electrode,
Foot switch,
2 m power cable
Technical data:
Maximum output power:
CUT cutting - 200 W-250 Ω,
BLEND mixed cutting - 150 W-200 Ω,
Forced coagulation FORCED COAG - 150 W-150 Ω,
Soft Coagulation SOFT COAG - 100 W-100 Ω,
BIPOLAR Bipolar Coagulation - 100 W-100 Ω,
Operating frequency - 360 kHz,
Patient electrical circuit - F,
Selection of supply voltage - 115-230 Vac,
Nominal frequency - 50-60 Hz,
Supply voltage - 350 VA,
Weight - 5.8 kg,
Dimensions in mm - 370x145x300,
Has a valid Technical Passport issued,
12 months warranty

