



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Diathermy Physioled THERMA (New)**

New

Italian production,

Therma is a medical device for temporary use, suitable for use in therapeutic diathermy capacitive procedures and resistive energy transfer of musculoskeletal injuries (contractures, sprains and strains, tendonitis, synovitis and bursitis) and osteoarticular pathologies, Therapeutic treatment of thermal and capacitive resistance is the latest and effective approach between diathermy treatments for all professionals who are looking for something more than a traditional diathermic device (microwave or walkie-talkies), Therma is a technological and therapeutic evolution that provides versatility, effectiveness and therapeutic benefits that are difficult to achieve with conventional diathermy devices. The treatment is similar to massage, during which you can provide the patient with energy that immediately feels a pleasant feeling of warmth,

Advantages:

A wide range of pre-set protocols,

Customized programs

Custom programs,

Capacitive mode for drainage and removal of contractions with a greater thermal effect,

Capacitive mode for drainage and relaxation of contractions with a fairly important thermal effect,

Resistive mode for biological stimulation with athermal effect or minimal heat,

Head with 3 working configurations

Technical data:

Maximum power: 150W / 250VA,

Operating frequencies: 500 kHz / 1 MHz,

Operating mode: Resistive / capacitive,

Power supply: 90-240 V / 50-60 Hz,

Electric input power: 300 VA,

Dimensions:

315mm wide,

Length 390 mm,

170 mm high,

Weight 5 kg,

Has a valid Technical Passport issued,

12 months warranty

