



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Physioled PRESSA compression massage apparatus (New)

New

Italian production,

Pressa is a modern and reliable new generation compression system. The equipment is intended for physiotherapeutic procedures in which mechanical pressure massage is desired for mechanical lymphatic drainage,

Accuracy of Press's measurement, multifunctionality of the planning system and operation of special physiological accessories, with overlapping sectors, allow to unblock lymph nodes before performing pressure massage,

The device was designed and manufactured using the latest electronic technology that supports up to 8 sectors,

Pressa devices use the automatic recognition system of the connected accessory, which avoids unwanted errors in choosing a protocol or program. Pressa allows to set different working pressure for each sector,

Advantages of the Physioled PRESSA compression massage apparatus:

Stored therapy protocols for each connected equipment,

Planning the sequence of treatment phases,

Very quiet compressor operation,

Automatic deflation at the end of the treatment,

Quick connection of accessories,

Technical data:

7 inch touch screen

Maximum output pressure 250 mmHg,

Pressure range in the sector 30-90 mmHg,

Number of sectors 8,

Power supply: 115 - 230 V

Electrical input 70 VA,

Treatment mode: Time or Cycles,

Dimensions:

315mm wide,

Length 390 mm,

170 mm high,

Weight 6 kg,

Has a valid Technical Passport issued,

12 months warranty

