



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Physioled LAMBDA 2 laser therapy apparatus (Monodiiodic + Exadiiodic) (New)

New,
Italian production,
LAMBDA 2 is laser therapy equipment,
Lambda is a programmable IR LASER control unit suitable for regenerative biostimulative and analgesic treatments,
It is possible to connect to optional lambda transducer terminals. Transmitters can have different dimensions and radiation power. This allows the operator to use the best handle according to the treatment chosen,
Lambda has a security control system that automatically deactivates the LASER emission if the emitting handle is bypassed or incorrectly positioned,
The Physioled LAMBDA 2 laser therapy apparatus has 2 laser heads 30W and 180W (6x30W) in the set
Advantages of Physioled LAMBDA 2:
A wide range of pre-set protocols,
Customized programs
Independent management of 2 different instruments ensured,
Head and six LASER emitters,
Automatic scanning at a certain frequency,
Visualization of supplied energy (Joule),
Warning signal on the head - it lights during power supply,
Automatic recognition of the connected accessory,
Technical data:
Output frequency: 0.5-10.000 Hz,
Head wavelength: 905 nm,
Output pulse length: 100 nsec,
Power supply: 90-240 V / 50-60 Hz,
Electrical input power (VA): 25 VA (lambda2),
Dimensions:
315mm wide,
Length 390 mm,
170 mm high,
Weight 3.5 kg,
Has a valid Technical Passport issued,
12 months warranty

