



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Physioled TESLA Plates Magnetotherapy device (New)

New

Italian production,

Tesla is a magnetotherapeutic ELF unit designed for professional applications,

Tesla produces a magnetic field (ELF) characterized by: low frequency, negligible electric field, low radiation power, negligible thermal effects, high tissue penetration,

Advantages of the Physioled Tesla Plates apparatus:

The camera has a graphic touch screen that allows interactive control of the device,

A wide range of pre-set protocols,

Customized programs

Automatic recognition of the connected accessory,

Phase sequence planning,

The device can be updated via USB

Technical data:

Output frequency: 2 - 100 Hz,

Maximum output level: 100/10 Gauss / mT,

Power supply: 90-240 V / 50-60 Hz,

Electric input power: 70 VA,

Dimensions:

315mm wide,

Length 390 mm,

170 mm high,

Weight 5 kg,

Has a valid Technical Passport issued,

12 months warranty

