



### **LIKO ROLLON patient lifter (Reconditioned)**

Used, very good technical condition,  
Swedish production (12-2014),

RollOn is a standing device designed for safe transfers to and from a sitting position. Perfect for situations requiring a transfer to or from the toilet, or between a bed and a chair. RollOn can also be an alternative to a wheelchair in some situations when the patient has an unbalanced balance, e.g. when he gets up for the first time after surgery or delivery. The RollOn stander works in such a way that the patient uses his strength to pull himself up to a standing position. When in the standing position, the patient is secured with protective and stable supports,

Simple, intuitive, manual and uncomplicated operation,

Safe and easy to open rear brackets,

Aluminum standing frame makes maneuvering easier,

Maximum load: 160 kg,

Height adjustable rear supports in the range of 96 - 114 cm,

Adjusting the distance between legs of the standing frame in the range of: 47.5 - 80 cm,

Wide range of height adjustment of front supports,

Wheels with a diameter of 75 mm (including two with brakes),

Dimensions: 62 x 92 x 140 cm,

Weight: 25 kg

User manual in English (PDF),

Has a valid Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month on the international market,

Financing option (Poland only): Installments, Leases, Loans,

