



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET 1002.21C0 side handle (Reconditioned)

Used, very good technical condition,
German production,
Used for lateral patient protection,
Better pressure distribution due to slightly concave shape,
Removable liner for easy cleaning,
Integrated clamping block,
Mattress dimensions: 9.5 x 21 x 4.5 cm,
Overall dimensions: 9.5 x 38 x 4.5 cm,
Has a valid Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month on the international market,
Financing option (Poland only): Installments, Leases, Loans

