



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SIEMENS AXIOM ARTIS DFC cardioangiograph (Reconditioned)

Used, very good technical condition,

German production,

The Siemens Axiom Artis platform is one of the most comprehensive interventional radiology systems on the market,

DFC model - floor mounted digital,

Arm C:

Rotation range:

LAO / RAO: 120 ° / 120 °,

CRA / CAU: 55 ° / 45 °,

Shoulder depth: 93 cm,

FDD range: 53 - 83 cm,

Range of arm support movement: from + 119 ° to -185 °,

Floor mounting,

Arm support weight: 665 kg,

Generator:

Nominal power: 100 kW,

High frequency generator,

Radiography: 50 - 125 kV,

Fluoroscopy: 50 - 125 kV,

Maximum lamp current: 1000 mA,

Fluorine pulse (frequency): 10, 15 and 30 pps,

Pulse width: from 3.2 to 10 ms,

Automatic fluoroscopy control (kV and mA),

High contrast fluoroscopy mode: 170 mA,

Type of MEGALIX Cat insert,

Focal dimensions: 0.4 / 0.8 mm,

Nominal anode input power: 35/80 kW,

Generator weight: 350 kg,

Table:

Longitudinal movement: 125 cm,

Lateral movement range: ± 17.5 cm,

Tilt within ± 15 ° (trendelenburg and lateral),

Height adjustment in the range: 77.5 - 110 cm,

Table top rotation range: ± 120 °,

Maximum patient weight: 200 kg,

Table top made of carbon fiber,

Top dimensions: 229 x 45 cm,

Weight: 452 kg,

Included:

Workstation - computers and monitors,

Software and full documentation,

Foot attachment,

Has a valid Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month on the international market,

Financing option (Poland only): Installments, Leases, Loans

