



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **The device for body shaping LPG ENDERMOLAB (Reconditioned)**

Used, perfect technical condition,  
French production (2018),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 830 W,

Endermology Alliance LPG

is a revolutionary technology that transforms the way we take care of our bodies. Thanks to the patented Endermologie Alliance, you can see significant results after just three treatments. This innovative method enables intense and rapid stimulation of tissues, while providing unparalleled comfort during the treatment,

The advanced technology affects the body in three dimensions, offering comprehensive care:

Fat reduction - Alliance endermologie effectively reduces localized and resistant fatty tissue, enabling body contouring without resorting to drastic methods,

Smoothing cellulite - fighting cellulite can often be a challenge, but with Alliance endermologie it becomes much easier. This advanced technology helps smooth the skin, reducing the appearance of cellulite and leaving it smooth and firm,

Skin firming - Alliance endermologie not only reduces fat and cellulite, but also stimulates the production of natural collagen, elastin and hyaluronic acid. As a result, the skin becomes more elastic, firm and youthful-looking,

Application:

Reducing secondary lymphedema of the arm after mastectomy,

Improving secondary lymphedema,

Improving lymphatic circulation in the treated area,

Relief of minor muscle pain,

Relief of muscle spasms,

Temporary improvement of local blood circulation,

Temporary relief of minor muscle pain associated with delayed muscle soreness,

Improving local circulation during burn rehabilitation,

Reducing the appearance of cellulite and reducing the circumference of treated areas,

Temporary improvement of lymphatic circulation and local blood circulation to improve skin trophicity in treated areas,

Improvement of skin quality, scarring, fibrosis,

Improvement of skin aging (wrinkles, fine lines, skin laxity, fat penetration, firmness, elasticity, complexion, bags under the eyes),

Stimulation of fibroblasts (synthesis of collagen, elastin, hyaluronic acid),

Radiation power:

Channel b: 12.46 dBm,

g channel: 9.33 dBm,

Channel n: 9.48 dBm,

Low-voltage 3.25 dBm,

Includes 4 heads:

Alliance 50 - intense and precise stimulation suited to any skin type, for smaller areas,

Alliance 80 - adapts to any tissue type, optimizes therapy intensity, for larger areas,

Motorized TR30 - for small and hard-to-reach treatment areas,

Ergolift - module for face and small scar type lesions, dynamically stimulates skin cells,

Treatment bed included,

Has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan,

