



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **USG SIEMENS ACUSON X300 (Reconditioned)**

Used, very good technical condition,

American production (2010),

Power supply: 230 V,

Frequency: 50/60 Hz,

The ACUSON X300 system is a portable imaging system designed for ultrasound diagnostics. The system uses advanced image processing technology,

The apparatus can be used in the following tests: abdominal, obstetric, gynecological, vaginal, urological, small organs (breast, thyroid testicles), musculoskeletal, orthopedic, cardiological, vascular tests,

The small size and mobility of the system allow access to the tested objects in conditions of limited space,

Imaging modes:

2D mode in double and split mode with THI basic and harmonic imaging,

M mode,

Doppler color mode,

Doppler power mode,

PW Doppler (Doppler pulse mode),

Equipped with:

Flat screen display,

Two active transducer ports and one rest port

DICOM function, possibility of connecting to the network,

Convex head - CH5-2,

Linear head - VF10-5,

Endovaginal head - EV9-4,

Illuminated control panel

Thermal printer SONY Model UP-D897,

USB port - for connecting an external printer,

RS-232C ports - for transferring data to a personal computer,

RGB, VGA output,

S-VHS in / out, Y / C in,

VHS In / Out (BNC),

Audio input / output

Accessory power socket

CD-RW / DVD-RW reader,

Transducer and gel holders,

Mobile on wheels with brakes,

Dimensions: 52 x 86 x 138/154 cm,

Cat 02,

Has a valid Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month on the international market,

Financing option (Poland only): Installments, Leases, Loans

