



### **Cosmetic chair with footrest (Reconditioned)**

Used, very good technical condition,

Swiss production,

230 V power supply,

50-60 Hz frequency,

110 W power

The negative pressure produced is -85 kPa / -638 mmHg,

Max capacity 30 l / min,

It has 1 MEDELA jar with a capacity of 2000 ml, (autoclavable) and MEDELA overflow jar with a capacity of 200 ml, (autoclavable),

Smooth suction regulation,

Built-in pressure gauge

Equipped with a strip for fixing the rings and additional equipment,

Includes a 2m long suction hose,

Mobile (on wheels with brakes).

Dimensions: 49 x 60 x 94 cm.

Cat 03,

Has a valid Technical Passport issued.

Warranty:

6 months for the domestic market - Poland,

1 month on the international market,

Financing option (Poland only): Installments, Leases, Loans

