



Enteral feeding pump SONDOMAT PLUS (Reconditioned)

Used, in very good condition,
German production,
Battery power supply - 4 x AA,
Flow rate: 1 - 300 ml/h,
Step of flow velocity: 1 ml/h,
Bolus speed: 1 - 500 ml/h,
Step of the bolus: 1 ml/h,
Delivery pressure: 175 kPa,
Battery alarm,
Dimensions: 8.7 x 16.6 x 4.5 cm,
Weight: 460 g,
He has issued a valid Technical Passport,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Instalments, Leasing, Loan

