



AESCLAP PROMIS LINE FLUID 2 laparoscopic pump (Reconditioned)

Used, in very good condition,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Electricity consumption: 70 W,

This is an advanced system that ensures maximum patient safety. Thanks to its simple and intuitive operation, it improves comfort and safety at work, and thus provides full functionality,

Maximum flushing and suction power: 500 mmHg,

Flow range: 0 - 1000 ml/min,

Dimensions: 305 x 135 x 328 mm,

Weight: approx. 10 kg,

Operating instructions in different languages (PDF),

He issued a valid technical passport,

Warranty:

6 months for the domestic market in Poland,

1 month for the international market,

Financing option (Poland only): Instalments, leasing, credit, loan

