



OLYMPUS CLV-FD light source (Reconditioned)

Used, in very good condition,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power consumption: 300 W,
Suitable for flexible endoscopes,
Xenon light - 300 W bulb output,
Emergency lighting - 150 W halogen bulb,
Air compressor with 3-stage adjustment,
Manual adjustment of light intensity,
Possibility of connecting a fiber optic cable with rigid optics,
Equipped with carrying handles,
Dimensions: 587 x 165 x 537 mm,
Weight: 22 kg,
He has issued a valid Technical Passport,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Instalments, Leasing, Loan,

