



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Arthroscopic pump OLYMPUS SRURGIPUMP (Reconditioned)**

Used, in very good condition,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power consumption: 70 W,

This is an advanced system that ensures maximum patient safety. Thanks to its simple and intuitive operation, it improves comfort and safety at work, and thus provides full functionality,

Maximum pumping pressure: 500 mmHg,

Flow range: 0 - 1600 ml / min,

Dimensions: 32.5 x 28.5 x 11.5 cm,

Weight: 7.3 kg,

He has issued a valid Technical Passport,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (only Poland): Instalments, Leasing, Loan

