



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRAISTECH diagnostic and treatment lamp (Reconditioned)

Used, in good condition,
Polish production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Bulb power: 60 W,
Flexible arm with length: 22 cm,
Bowl diameter: 20 cm,
Manual height adjustment in range: 49 cm,
Maximum height: 200 cm,
Base diameter: 30 cm,
He has issued a valid Technical Passport,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (only Poland): Instalments, Leasing, Loan

