



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Pulse oximeter SMITHS MEDICAL BCI 3301 (Reconditioned)**

Used, in good working order,  
American production,  
Battery power supply: 3 x R14,  
Durable and lightweight construction,  
Ergonomically designed to fit comfortably in your hand,  
Bright and easy to read display,  
Quantitative pulse force bar on the screen,  
Fast, reliable measurements of SpO2 and heart rate,  
Stores data for up to 99 different patients,  
Simple on/off control with automatic off to save battery life,  
Designed for patients of all ages, from newborns to adults,  
Included:  
Clip sensor with extension cable,  
Batteries 3 x R14,  
Dimensions: 8.5 x 17.5 x 4 cm,  
Weight: 0.4 kg,  
He has issued a valid Technical Passport,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (only Poland): Instalments, Leasing, Loan

