



### **Pulse oximeter SCHILLER OXIMETER OX-2 with docking station (Reconditioned)**

Used, very good technical condition,  
U.S.A. production.

Power supply: battery (replacement at extra charge), and through docking station - 230 V,

Frequency: 50/60 Hz,

SCHILLER OX - 2 provides the most precise and reliable SpO2 readings, also during movement and in states of decreased perfusion,

Pulse oximeter with detachable panel/screen for patient transport allowing 4 hours of operation thanks to built-in rechargeable battery,

Masimo SET signal extraction technology provides excellent measurement quality in patient motion and low perfusion conditions,

Provides SpO2 and heart rate measurement:

SpO2 measurement range: from 0% to 100%,

Pulse measurement range: from 25 bpm (beats per minute) to 240 bpm,

Perfusion range: from 0.02% to 20%,

Built-in display,

Audible and visual alarms,

Selectable menu language: EN, DE, IT, SL, FR, NL, DA, PT, ES,

Volume control and MUTE function,

Automatic calibration of the device,

Included:

Docking station,

Masimo clip sensor,

Power cord,

Dimensions: 27.5 x 9 x 19 cm,

Weight: 1.8 kg,

Has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for international market,

Financing possibility (Poland only): Installments, Leasing, Loan

