



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **LEF welding machine - LABO ELECRTOFRANCE TS47TUSB (Reconditioned)**

---

Used, in very good condition,  
French production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power consumption: 500 W,  
Very quick welding,  
Barcode editing software,  
Printer with high quality 1 line printout,  
Weld width: 15 mm,  
Welding speed: 10 m/min,  
Set the temperature in the range: 20 - 230°C with an accuracy of 1%,  
Colour touch screen,  
Intuitive multifunctional menu,  
Counter of reset weld seams,  
Paper guide - adjustable depth in range: 0 - 30 mm (from the edge of the bag),  
Automatic start as soon as the set temperature is reached,  
Made of stainless steel,  
USB output,  
Included:  
Barcode reader,  
USB stick with software,  
Roller conveyor,  
Dimensions of the welding machine: 550 x 260 x 180 mm,  
Dimensions of the roller conveyor: 800 x 300 x 60 mm,  
Weight of the welder: 12 kg,  
He has issued a valid Technical Passport,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (only Poland): Instalments, Leasing, Loan

