



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Vacuum generator EUROSETS MOBI.S (Reconditioned)**

Used, in good working order,  
Italian Production, (2006),  
Power supply: 230 V or battery (2 x AA),  
Continuous electronic control via microprocessor,  
The current pressure and setpoint are always shown in the display,  
Time controller,  
3 modes of operation:  
Aspiration by oscillation in range,  
Continuous Aspiration ON-OFF,  
Aspiration by value,  
Selection of the displayed pressure unit: CmH<sub>2</sub>O, mmHg, mBar, Pax100,  
Acoustic and visual alarms,  
Menu languages: Italian, English, German,  
DMS - data management system,  
Included:  
Original suitcase,  
Power supply,  
Cat 02,  
He has issued a valid Technical Passport,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (only Poland): Instalments, Leasing, Loan

