



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment table Wesseling KOMBI ELITE S Hydraulic (New)

New guy,
Dutch production,
The Elite S Combo is the sister of the successful Elite Combo Combo.
This table has a shorter part of the seat, 40 cm long, so that the patients are better prepared for the treatments.
This table has a large height adjustment range for an ergonomic working environment.
Due to the highest position of 98cm, treatments can often be performed in a standing position.
Due to the unique separation of the sections, this table can be used for almost any treatment in general or gynaecological practice.
Round legs and thicker table profiles not only look beautiful, but are also very strong and stable.
The table height can be adjusted by means of a hydraulic cylinder,
Technical parameters:
Length: 188 cm,
Width: 64 cm (Optional 70 cm or 80 cm)
Number of sections: 2,
Height adjustment: from 48 to 98 cm,
Adjustable headrest angle: $-35^{\circ}/+85^{\circ}$
Back section adjustment: $0^{\circ}/+85^{\circ}$.
Wheels: 75mm diameter
Maximum load: 200kg,
Additional options are available (at an extra charge):
Grounding point,
Memoryfoam,
Special width of tabletop (70cm),
Special width of tabletop (80cm),
Additional footswitch,
Adapters for shoulders/feet/footrests,
Crotch incision,
Silicium (set),
A set of footrests,
Footrests,
Earthing cable,
Table height adjustment system
Paper roll holder,
Headrest pillow,
Plug for two switches,
Nose opening,
Armrest
Additional 10 000N actuator with 300mm stroke,
Additional 12 000N actuator with 300mm stroke,
Additional high-speed 8 000N actuator with 300mm stroke,
Thicker upholstery (70mm),
Additional manual switch,
He has issued a valid Technical Passport,
Warranty 24 months,
Financing possibility (Poland only): Instalments, Leasing, Loan

