



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **Pressure relief mattress with SELECT AIR MAX PUMP 240 variable pressure pump (Reconditioned)**

---

Used, in very good condition,  
of the U.S.A. Production Department,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 720 W,  
Quiet, lightweight, economical and reliable system,  
Helps to manage the skin's microclimate and improve moisture evacuation,  
The hard mattress mode facilitates nursing care during dressing changes and patient repositioning,  
Vertical mode provides additional support when the upper part of the bed is raised,  
10 levels of mattress hardness,  
CPR quick exhaust function for emergency procedures,  
3 levels of mattress pulsation power,  
The foot cushion helps to cushion the foot,  
The zippered cover is impervious to fluids and vapour-permeable,  
Individual partitions divided into 3 sections: head, seat and foot,  
Integrated side airbags with foot cushion,  
Maximum load: approx. 180 kg,  
Dimensions of the pump: 35.6 x 17.8 x 15.2 cm,  
Weight of the pump: 4.1 kg,  
Mattress weight: 3.2 kg,  
Operating instructions in Swedish,  
He has issued a valid Technical Passport,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (only Poland): Instalments, Leasing, Loan

