



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Coagulology analyzer HELENA C-1 (Reconditioned)**

Used, in very good condition,  
English production,  
Power supply: via 12 V power supply,  
Frequency: 50/60 Hz,  
Power consumption: 12 W,  
Semi-automatic,

The innovative line of C-series coagulometers combined with a wide panel of Helena reagents ensures high quality results,  
Multiple measurement options - autostart sensors, minimum reagent and sample volumes, optional software for performance analysis and data management,

A complete, simple and economical solution for any laboratory,  
All C-series coagulometers perform coagulation, chromogenic and D-Dimer tests,  
One measurement channel,

6 thermostatic sample positions,

2 thermostatic reagent positions,

A reaction vessel in the form of a single cuvette,

Double-sided markings and average values,

Results shown on the display,

3-point calibration curve,

Reagent consumption: 0.25 to 0.5 µl,

Included:

12 V power supply unit,

User manual in Polish,

Dimensions: 245 x 130 x 60 cm,

Weight: 0.51 kg,

He has issued a valid Technical Passport,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (only Poland): Instalments, Leasing, Loan

