



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Spirometer MIR SPIROLAB (Reconditioned)

Used, in very good condition,
Italian production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power consumption: 20 VA,
No calibration required,
The device measures the following parameters: VC, ERV, FEF25, FEF50, FEF75, FEF25/75, FET, FEV1, FEV1%, FIV1, FIVC, FIV1%, FVC, IVC, MVVV, PEF, PIF, EXT,
Easy to read and easy to read control panel,
Colour LCD display: 16 lines x 20 characters (128 x 160 pixels),
Integrated thermal printer with a paper width of 110 mm,
Printing width: 80 mm,
Possibility of connecting to a PC via RS232 interface,
Flow rate: +/- 20 l/s,
Flow accuracy: +/- 5% or 200 ml,
A set is also included in the set:
Control panel,
Transmitter with cable,
Adapter for turbines,
Reusable turbine,
Power cord,
User manual in English (PDF),
Dimensions: 350 x 280 x 112 mm,
Weight: 4.9 kg,
Cat 01,
He issued a valid technical passport,
Warranty:
6 months for the domestic market in Poland,
1 month for the international market,
Financing possibility (Poland only): Instalments, Leasing, Loan

