



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Light source AOHUA LG150-2 (Reconditioned)

Used, good technical condition,

P.R.C. Production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 250 VA,

It has local scratches of paint on the housing,

Design with double halogen bulb 15 V / 150 W,

Switching the light source on the front panel,

Ten-step intensity control,

Illumination: ≥ 70000 Lx,

Colour temperature: ≥ 3200 K,

Air pump pressure: 30 - 60 kPa,

Airflow: 4 - 10 l/min,

Integrated air water pump for gastroscopy,

Water bottle holder,

Power cord included,

Dimensions: 38 x 41 x 14.5 cm,

Weight: 11 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

