



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED shank holders WS-04.5 (Reconditioned)

Used, very good condition,

Polish production,

They serve to immobilize the patient's lower limbs in a lying position,

The single handle consists of a polyurethane mattress which is a shank clamp and a support part made of stainless steel for height adjustment,

Can be retrofitted with a rail mounting clamp at additional cost,

Working length of the support part: about 40 cm,

Mattress dimensions: 19 x 6.5 x 3 cm,

Handle weight: 1.6 kg,

The weight of the set: 3.2 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan,

