



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SUNRISE MEDICAL STERLING SWIFT (Ex-Display)

Used, very good condition,
Power supply: rechargeable - two new batteries included,
Easy to use and very agile 3-wheel electric scooter,
Thanks to its compact dimensions, it is suitable for all road conditions,
He is doing well in urban conditions as well as on field, forest roads,
Very strong brakes ensure safety,
Fully amortized suspension,
Rear tyres trained,
Large, comfortable, adjustable armchair with folding armrests,
Removable basket on the steering wheel,
Height-adjustable ball holder,
Full road lighting:
Position lights,
Indicators,
Stop light,
Maximum speed: 12 km/h,
Range: up to 30 km (33 Ah batteries according to manufacturer),
Maximum patient weight: 136 kg,
Wheel diameter: 27 cm,
Clearance: 10 cm,
Dimensions:
Seat: 51 x 45 x 11 cm,
Backrest: 46 x 38 x 8 cm,
Total width: 61 cm,
Total length: 132 cm,
Total height: 103 cm,
Battery charger included,
He has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Instalments, Leasing, Loan

