



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LID Lisa LED light (New)

- New,
- Made in France,
- LED reading light with wired remote switch, popular for use in patients' rooms. A second, easy-to-use switch (low safety voltage 12 V) connected to a 150 cm long cable allows bedridden or mobility-impaired people to use the lamp easily and intuitively during the day and night. It can also be used as a bedside lamp.
- It provides soft, uniform, natural light (40° beam width) that does not cause glare or micro-flicker that can cause headaches and eye fatigue,
- The sturdy 60 cm long gooseneck makes it easy to adjust the lamp head even with one hand,
- Specifications:
 - Light intensity: 420 Lumens,
 - Light color: 4500°K* (neutral white),
 - Color rendering index: IRC 95/100,
 - Electricity consumption: 3 W,
 - LED lifetime: minimum 50,000 hours,
 - Power at a distance:
 - 50 cm: 3500 lux*.
 - 30 cm: 10400 lux*.
 - Electrical protection: insulation class II,
 - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
 - Power cord length: 190 cm,
 - Available colors : white,
- Dimensions:
 - Lamp head dimension: 60 x 100 mm,
 - Flexible arm: 600 mm,
 - Weight: 1.2 / 2.6 kg (with metal base),
- Possible attachment: to the base or clamp,
- Has a current Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

