



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LID Thera DUO Infrared light (New)

- New,
- Made in France,
- Dual-head infrared lamp with two 250 W bulbs, on two 50 cm long goose necks, which are very durable and quiet. A session can be performed in half the time on two different zones (neck/arms - neck/back - both thighs, etc.) with better heat distribution than with one lamp, or one treatment zone can be treated more intensively (500 W) with careful placement of both heads,
- The light can also be used to affect all sides of a muscle group (e.g. biceps femoris, rectus femoris, vastus lateralis...). Using only one plug, Thera Duo is functionally mobile and occupies a small space,
- Recommended for specialties such as physiotherapy / massage / acupuncture / orthopedics,
- Robust design allows for accurate and reliable positioning,
- Ergonomic handle allows precise positioning,
- Specifications:
 - Electricity consumption: 250 W x 2,
 - Wavelength: 1200 nm,
 - Class 1 medical device,
 - IR bulb life (E27): 5000 hours minimum,
 - Electrical protection: insulation class I,
 - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
 - Power cord length: 4 meters,
 - Available colors : white,
- Dimensions:
 - Dimension of the lamp head : 260 mm diameter,
 - Dimension of the base: 69 x 69 x 115 mm,
 - Weight : 4 kg,
 - Length : Flexible arm: 2 x 50 cm,
- For an extra charge, you can retrofit the lamp with a timer switch,
- Additional mounting accessories are required for the lamp , of your choice (at an extra charge):
 - Telescopic stand base 60-100 cm white metal, heavy,
 - Telescopic stand base 60-100 cm white metal, heavy -version with brakes in 2 wheels,
 - Telescopic stand base 90-130 cm metal white, heavy,
 - Telescopic stand base 90-130 cm metal white, heavy -version with brakes in 2 wheels,
 - Plain clamp - black,
 - Clamp for medical rails with vertical lamp outlet,
 - Universal clamp,
 - Wall bracket,
 - ARM wall bracket white (2x40cm),
 - ARM wall mount white (1x40cm),
 - Tabletop bracket ARM white (1x60cm),
 - Plate with metal pipe in the middle $\varnothing 16$,
 - Plate with metal tube on the edge $\varnothing 16$,
 - Plate on ARM holder,
 - Female socket with adaptor plate $\varnothing 16$ mm x 27 mm,
 - Adapter plate male rod $\varnothing 15$ mm x 30 mm,
- It has a valid Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

