



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Telescopic stand base - heavy duty (New)

- New,
- Made in France,
- Base in the shape of a five-pointed star with a diameter of 53 cm with wheels (diameter 5 cm). Very stable,
- Frame made of easy to clean and disinfect smooth steel coated with epoxy resin. Castor base without indentations or grooves, facilitating effective cleaning and reducing the retention of dust or other debris (ideal for hospitals and health centers),
- Adjustable column height, allowing the lamp to be adjusted to the appropriate height for all medical applications,
- Standard stand height is 60-100 cm, for an extra charge you can purchase a 90-130 cm stand,
- For an extra charge, it is possible to purchase a stand with additional brakes in 2 casters,
- Dedicated to Bella, Bella Lumax, Thera and Thera Duo lights,
- Height: 60 - 100 cm / 90 - 130 cm,
- Weight: 5.5 / 5.8 kg,
- Available colors : white or black,
- Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

